

Cryotainer B.V.
Fokkerstraat 507
3125 BD Schiedam
Phone: +31 (0)10 262 30 53

P.O. Box 12086
3004 GB Rotterdam
www.cryotainer.com



Datasheet FT 1500 (Flying Tainer)



Dimensions l x h x b	mm	12.200 x 3.650 x 2.500
Maximum working pressure	bar	22
Maximum flow	m ³ /h	1.500 (2 x 750)
Maximum outlet temperature nitrogen	°C	15 ⁰ till 20 ⁰ below ambient temperature
Total weight	kg	38.000
Suitable for cryogenic gasses		Liquid Nitrogen; Liquid Oxygen*
Storage capacity	ltr	19.700
Outlet connection		Both sides 2" weco fig. 100

* For use of Liquid Oxygen Container and Vaporiser should be purged for use of it