

Cryotainer B.V.  
Fokkerstraat 507  
3125 BD Schiedam  
Phone: +31 (0)10 262 30 53

P.O. Box 12086  
3004 GB Rotterdam  
www.cryotainer.com



## Datasheet FT 1700 (Flying Tainer)



Dimensions l x h x b	mm	12.200 x 3.650 x 2.500
Maximum working pressure	bar	22
Maximum flow	m <sup>3</sup> /h	1.700 ( 2 x 850 )
Maximum outlet temperature nitrogen	°C	15 <sup>o</sup> till 20 <sup>o</sup> below ambient temperature
Total weight	kg	38.000
Suitable for cryogenic gasses		Liquid Nitrogen; Liquid Oxygen*
Storage capacity	ltr	19.700
Outlet connection		Both sides 2" weco fig. 100

\* For use of Liquid Oxygen Container and Vaporiser should be purged for use of it

Chamber of Comm. no. 24361740  
VAT no. NL 813168879 B01  
Bank: Fortisbank 97.10.52.301  
IBAN: NL28.FTSB.0971.0523.01  
BIC: FTSBNL2R